

WORKSHOP PRESENTATION

Open Access

High-resolution three-dimensional ANGIE T1 mapping of the whole heart

Bhairav B Mehta^{1*}, Michael Salerno^{2,1}, Frederick H Epstein^{1,3}

From 18th Annual SCMR Scientific Sessions Nice, France. 4-7 February 2015

Background

Assessment of fibrosis in thinner myocardial structures such as the right ventricular (RV) and left atrial walls would be valuable in disorders such as pulmonary hypertension, congenital heart disease, and atrial fibrillation. We recently developed a novel technique, termed ANGIE, which provides high-resolution 2D T1 mapping for assessment of thin structures such as the wall of the RV. However, high resolution 2D imaging is limited by lower SNR and time efficiency. A 3D image acquisition provides higher SNR and complete coverage of the heart. However standard 3D navigated techniques have prohibitively long scan times. Thus, the aim of the present study was to extend ANGIE to perform high-resolution three-dimensional (3D) T1 mapping of the whole heart within a clinically acceptable scan time.

Methods


A novel image reconstruction algorithm, which uses compressed sensing (CS) with local low-rank sparsity in conjunction with partial Fourier and parallel imaging, was developed to achieve high acceleration (R=10) for 3D ANGIE. High-resolution (1.4×1.4×4mm³) T1 mapping of the whole heart was performed using 3D ANGIE in three healthy volunteers on a 1.5T system (Avanto Siemens). Imaging parameters included: matrix size=224×224×26, number of inversion times=12, and partial Fourier reduction factor of 3/4 along phase and partition encode. Additionally, 2D MOLLI (2.5×2×4mm³) and 2D ANGIE (1.3×1.3×4mm³) scans were performed at three different short-axis slice positions.

Results

Figure 1 illustrates example reconstructed (A-D), fully sampled (E-H) and reconstruction error (I-L) images

 $\overline{}^{1}\text{Department}$ of Biomedical Engineering, University of Virginia, Charlottesville, VA, USA

Full list of author information is available at the end of the article



Conclusions

3D ANGIE provides accurate high-resolution $(1.4\times1.4\times4\text{mm}^3)$ native T1 maps of the whole heart within a clinically acceptable scan time.

Funding

NIH R01 EB 001763, AHA Grant-in-Aid 12GRNT12050301, AHA 14PRE20210008, and Siemens Medical Solutions.

Authors' details

¹Department of Biomedical Engineering, University of Virginia, Charlottesville, VA, USA. ²Medicine (Cardiology), University of Virginia, Charlottesville, VA, USA. ³Department of Radiology and Medical Imaging, University of Virginia, Charlottesville, VA, USA.






Figure 1 Example reconstruction results from a retrospectively undersampled (rate 10) fully sampled 3D ANGIE dataset. Reconstructed (A-D), fully sampled (E-H) and corresponding reconstruction error (I-L) images from a single partition at four inversion time points of a 3D ANGIE dataset acquired from a healthy volunteer. The image reconstruction algorithm effectively suppresses aliasing artifacts and retains image spatial resolution. These results illustrate the capability of achieving high acceleration rate (R=10) using the proposed image reconstruction algorithm.




Figure 2 Comparison of T1 mapping techniques in healthy volunteer. Example native T1 maps from a healthy volunteer at three different slice positions acquired using 3D ANGIE (A-C), 2D ANGIE (D-F) and 2D MOLLI (G-I). 2D MOLLI provides high SNR but incomplete definition of the RV wall, whereas 2D ANGIE provides good definition of the RV wall but lower SNR. However, 3D ANGIE provides both good definition of the RV wall as high SNR.

Mehta et al. Journal of Cardiovascular Magnetic Resonance 2015, **17**(Suppl 1):W4 http://www.jcmr-online.com/content/17/S1/W4

Published: 3 February 2015

doi:10.1186/1532-429X-17-S1-W4

Cite this article as: Mehta et al.: High-resolution three-dimensional

ANGIE T1 mapping of the whole heart. Journal of Cardiovascular Magnetic

Resonance 2015 17(Suppl 1):W4.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit

